

# CheckCLIP

Single-Use | Medication ID



## A *simply safer* way of labelling syringes.

CheckCLIP safely attaches the ampoule to the syringe for accurate cross check without obscuring incremental dosage marks on the syringe to reduce medication errors.



*Simply clip together*

# Mistake-proof the labelling of syringes.

Quicksmart CheckCLIP is a simply safer way to label syringes with the name of the medication that has been drawn up always visible. This significantly reduces the risk of morbidity and mortality caused by medication errors.

## Patient safety

Traditional labelling techniques that involve handwriting on the syringe, attaching the ampoule with surgical tape or labelling with stickers are all prone to medication errors. Tape and stickers obscure incremental marks on the syringe and handwriting is time-consuming and open to transcription error. As a result, virtually every practitioner experiences a 'near miss' or worse at some stage in their chaotic work environment.

CheckCLIP avoids these traps by safely attaching the ampoule to the syringe for accurate cross checking without obscuring incremental dosage marks on the syringe. It also covers any exposed sharp ampoule edges. With Quicksmart CheckCLIP the ampoule containing all the drug details remains with the syringe at all times.

## Staff safety - physical and psychological

Medication error poses serious risks to patient safety, regulatory compliance and hospital litigation exposure. It is also a significant occupational health and safety concern for practitioners. One third of all doctors who have experienced a 'near miss' medication error report that the stress and repercussions negatively impacted on their life. 97% of nurses worry about medication errors. **68% of health professionals believe medication error can be reduced by better labelling. CheckCLIP provides that better labelling solution.**

A study by the Joint Commission of 3,171 physicians who had experienced (committed) a medication error found that<sup>6</sup>:

- 61% reported increased anxiety about future errors
- 44% reported a loss of confidence
- 42% reported sleeping difficulties
- 42% reported reduced job satisfaction
- 13% reported harm to their reputations


## Medication and syringe swap errors

Medication errors kill at least 1 patient every day in American hospitals and another 1.3 million people are injured every year at a cost of \$6 billion annually.<sup>1,2</sup>

A study conducted by Department of Anesthesia, Sunnybrook and Women's College Health Science Centre, University of Toronto, Ontario, Canada showed that 70.4% of 687 anaesthetists surveyed indicated syringe swap is one of the most common contributing factors to medication errors.<sup>3</sup>

A South African study found 94% of participants (125 of 133) admitted to having inadvertently administered a wrong drug. Contributory causes identified included<sup>4</sup>:

- Syringe swap - 40%
- Misidentification of drugs - 27.1%
- Fatigue - 14.1%
- Distractions - 4.7%
- Mislabelling of syringes - 4.7%




## Common factors in medication errors

Study of Injectable Medication Errors done by the American Nurse Association and Inviro Medical in 2007 showed the following results<sup>6</sup>:

- Too rushed / busy environment - 78%
- Poor / illegible handwriting - 68%
- Missed or mistaken physicians orders - 62%
- Similar drug names or medication appearance - 56%
- Working with too many medications - 60%
- Only 37% claim Injectable Medications are always labelled

Factors Contributing to Medication Errors<sup>4</sup>



1. Harden, S. (2006). "One Death Per Day Due To Medication Errors: Hospitals Seek Improvement By Learning From The Best In Aviation ". Press Release Newswire.  
 2. Kohn, L. T., J. Corrigan, et al. (1999). To Err is Human - Building a Safer Health System, IOM HealthGrades Study, National Academy Press.  
 3. Orser, B. A., R. J. Chen, et al. (2001). "Medication errors in anesthetic practice: a survey of 687 practitioners." Can J Anaesth 48(2): 139-46.  
 4. Gordon, P. C., R. L. Llewellyn, et al. (2006). "Drug administration errors by South African anaesthetists--a survey." S Afr Med J 96(7): 630-2.  
 5. Waterman AD, G. J., Hazel E, et al (2007). "The emotional impact of medical errors on practicing physicians in the United States and Canada." Jt Comm J Qual Patient Saf 33: 467-476.  
 6. "2007 Study of Injectable Medication Errors - An Independent Study Sponsored by the American Nurses Association and Inviro Medical Devices." from <http://www.nursingworld.org/MainMenuCategories/OccupationalandEnvironmental/occupationalhealth/SafeNeedles/2007InviroStudy.aspx>

## Simply attach the ampoule to the syringe

CheckCLIP enables the ampoule to travel with the syringe for visual cross-checking and accurate administration to the patient. Simply:

- Enables visual cross checking of drug and dose details
- Ensures a clear unobscured view of syringe increments
- Prevents sharps injuries from opened glass ampoules
- Provides clear indication of drug expiry dates
- Avoids accidental ampoule swaps and syringe swaps



## Mistake-proof process

Without ever letting go of the ampoule or syringe, you can:

- ATTACH the disposable CheckCLIP to the syringe
- SNAP the ampoule open and draw up the drug
- CLIP the ampoule safely to the CheckCLIP
- CHECK drug and dose with a clear view


There is no chance of syringe or ampoule swap error. By merely attaching then clipping, you not only save time but also lives in emergency and routine situations.

**Remember:** *double check and be sure.*

## Applications of CheckCLIP

CheckCLIP should be used wherever drugs are drawn into a syringe from an ampoule or vial.

- Hospitals - including Emergency Department, Operating Room, Anaesthetics
- Paramedics
- Medical Centres
- Animal labs



CheckCLIP vs. Current Practices	CheckCLIP	Medical Sticky Tape	Hand Written Label	Drug Class Colour Coded Label
Attach ampoule to syringe	Yes	Yes	No	No
Visual cross checking of drug/dose details	Yes	Poor or obscured	No	No
View of syringe increments	Yes	No	No	No
Clear indication of drug concentrations	No	No	Yes, if written on label correctly	No
Prevent sharps injuries from glass ampoules	Yes	No	No	No
Clear indication of drug expiry date	Yes	No, obscured by taping	No	No
Risk of mix-up > Ampoule leaves hand during process	<b>Almost nil &gt; Ampoule does not leave hand</b>	Potential > Ampoule leaves hand	Potential > Ampoule leaves hand	Potential > Ampoule leaves hand

# Specifications


Product edition	For Ampoule	For Vial
Product code	QSCCGEN	QSCCVIAL
Compatible syringe sizes	1-20ml luer and luer lock (majority)*	
Compatible ampoule / vial sizes	Glass and plastic ampoules 1-10ml (majority)**	Vials with diameter of up to 28mm*** (The sum of the syringe radius and vial radius must be less than 23mm)
Size (per unit)	Length 3.7cm Width 2cm Height 2.5cm	Length 4.2cm Width 2cm Height 1.8cm
Weight (per unit)	Approximately 1g	Approximately 1g
GMDN code	46240	

\* Compatibility depends on geometry and diameter of the syringe  
 \*\* Exceptions include ampoules 25ml or larger  
 \*\*\* Compatibility depends on the geometry and diameter of the vial


Material	Borealis polypropylene (PP)
Latex	No
DEHP	No
Sterile	No
Re-usable	No. Each unit attaches 1 ampoule or vial onto a syringe.
Disposal method	Once medication has been administered and drug identification is no longer required, the whole syringe unit (with CheckCLIP and ampoule or vial attached) is disposed of as medical waste in a sharps container.
Regulatory approvals	Australia - Listed under TGAARTG NO: 134522 USA - FDA listed

## Instructions for Use


### CheckCLIP FOR AMPOULE




1. Face barb down and push CheckCLIP on to syringe using your thumb.



2. Snap open the ampoule and draw up the prescribed medication. Without the ampoule or syringe leaving your hands, proceed to step 3.



3. Push ampoule up on to CheckCLIP.





4. Double check before injecting.


## CheckCLIP

Single-Use | Medication ID


### CheckCLIP FOR VIAL




1. Face barb down and push CheckCLIP on to syringe using your thumb.



2. Draw up the prescribed medication and without the vial or syringe leaving your hands, proceed to step 3.



3. Using the needle track as a guide, push vial up on to CheckCLIP.



4. Double check before injecting.

**QLICKSMART Pty Ltd.**  
 +61-7-3844-1182 T  
 +61-7-3844-1183 F  
 mjs@quicksmart.com  
 www.quicksmart.com

