



Work Safe

Play Safe Tomorrow.

▲ The underestimated cut injury from glass ampoule

One in three healthcare workers has experienced a cut when opening a glass ampoule by hand. Regrettably, many users and clinical administrators have underestimated the likelihood of these injuries.

Cuts from glass ampoules can happen to anyone, including you, no matter how experienced you are in opening ampoules. Even a minor ampoule cut can restrict your clinical activities. More serious injuries such as tendon lacerations have been reported.

Repeated personal injuries whilst opening ampoules, led an Australian registered nurse to invent the award-winning Qlicksmart SnapIT ampoule opener.

▲ The award-winning ampoule opener



SnapIT is designed to protect users from injury when an ampoule exposes a sharp edge or shatters while being opened.



SnapIT ampoule openers are made from sturdy, high-quality plastic.



Being reusable, SnapIT is a cost-effective solution to open glass ampoules.




SnapIT is designed to allow for customization with your institution's name and logo, which is a great way of reinforcing to your customers and staff that their safety is important to you.



The ampoule openers are simple to use with a quick "insert and snap" action.



The ampoule opener has won multiple awards, which demonstrate its practicality and effectiveness in reducing glass ampoule injuries.





SnapIT Lite is perfect for hanging on a lanyard or carrying in a pocket for easy, on-the-go use.

Workplace safety is the right of all healthcare workers, and the responsibility of all employers.

Use the Qlicksmart SnaplT ampoule openers and be safe.

■ Specifications



Product edition	Lite	
Product code	SN-02R	SN-02L
Product size	Regular	Large
Compatible ampoule sizes*	1-2ml, 5-10ml, 10-15ml	5-10ml, 10-15ml, 20-25ml
Size per unit	Length: 80mm Width: 13mm	Length: 80mm Width: 18mm
Weight per unit	11g	13g
Colour range	Metallic Grey, Purple, Teal, Pink, Blue.	Burgundy
GMDN code	10098	
Material	Plastic (body), Silicon (O-ring), Stainless steel (keyring and spring).	
Latex	No	
DEHP	No	
Sterile	No, not autoclavable.	
Re-usable	Multi-use. O-ring can be replaced if lost or broken.	
Regulatory	USA FDA listed	


*Compatibility depends on the ampoule's geometry and the diameter of the ampoule lid.

Ampoule of the same size can have a different geometry and amopoule lid diameter. The ampoule size in ml is only used as a guideline.


Snapit

Multi-Use Ampoule Opener

▲ How to use SnapIT Lite




 Align the dots if present. Gently insert all of the ampoule lit into the Snapit Lite. The ampoule lid becomes locked in at the neck.



2. With little effort, lever away from the dot to snap off the lid. The lid should remain held in the Snapır.



Carefully aim over a sharps bin and eject the lid by pushing down on the rod with your thumb.



4. Clean your Snapit regularly. Remove O-ring as shown to disassemble the parts and clean with water or hospital grade cleaner. Do this in a clean dish to avoid loss of parts. Dry completely then re-assemble.

NOTE: If the ampoule lid breaks, tap the Snapit Lite over a sharps bin to remove excess glass pieces then rinse under running water before disassembling.

SIZES

If the ampoule lid is too big, it will not fit through the entrance. Use a larger Snapı τ . If the ampoule is too small, the Snapı τ will not hold the top securely. Use a smaller Snapı τ .



Qlicksmart Pty. Ltd.



Distributed by:



Level 1, 148 Boundary Street, West End, Brisbane, QLD 4101, Australia



+61-7-3844-1183



www.qlicksmart.com

+61-7-3844-1182



hello@qlicksmart.com

Stay Connected









